


AMRAD  **Medical**



AAU ELITE


U-ARM RADIOGRAPHIC AND INTEGRATED DR SYSTEM



Large display of arm and detector rotation



Large LCD display with easy touch buttons for accurate positioning



Near-patient controls offer efficient equipment positioning refinements



Position equipment from anywhere in the room



Choose from a variety of mobile tables to add to your system. Available, if needed.

Responsible Dose Management

Low Dose Sensitivity Mode (LDSM) allows for up to 50% dose reduction in dose sensitive studies while maintaining excellent diagnostic quality.

Peak Performance

Achieve peak performance by migrating most exams to standing positions improving productivity and throughput.

Seamless Stitching

Perform LLI and Scoliosis exams with accuracy and ease using the automated equipment and automatic stitching function.

Smart Investment


AAU Elite offers equal features and earn equal reimbursement to all of its more costly competition and all of this at 5 times the warranty. Smart decisions equal smart investments. Be smart!



Refine equipment positions and remove complicated equipment angles with conveniently located control panel

**Improve Workflow
with protocol-driven
movements.**

AMRAD® Medical DR Elite™ Image Capture Workstations



- Simple, clean, and intuitive two screen workflow: (Work list and Acquire/Review screen)
- Multi-touch or keyboard interface with touch gesture and animated sliding transparencies
- Image preview in 2 seconds from exposure
- Integration with Electronic Medical Record Systems (EMS) with full MWL, MPPS integrations with RIS/HIS providers
- Multiple automated exports that can be user configured with a single touch of a button
- Complete integration of dose area product meters (DAP) with all readings recorded for export as well as embedded in DICOM image tags

Advanced Features

- Generator Interface Control
- Auto-Positioning
- Low Dose Sensitivity Mode (LDSM) allows for as much as 50% dose reduction in dose sensitive studies while maintaining excellent diagnostic image results
- Rejection analysis with user customization of reports
- Embedded Operating System for network safe environments
- Embedded Help Functions for immediate support services
- Windows 10 Embedded Operating System

Options

- Stitching
- Workstations include 24/7/365 Support and 99% Uptime Guarantee



www.amradmedical.com

SUMMIT INDUSTRIES, LLC
7555 N. Caldwell Avenue, Niles, Illinois 60714
(p) 773-353-4000 | (p) 1-800-729-9729



LITR041-1120

SAI GLOBAL
ISO 13485
Quality