

AMRAD Medical




CS Classic

WALL MOUNTED CHEST STAND SYSTEM

The AMRAD Medical CS system is designed for facilities that specialize in COVID-19 screening, Pulmonology and Occupational Medicine.

- 30-50kW generators
- Extensive vertical travel – 14”-73”
- Accommodates a full range of upright exams
- Pair with optional mobile table to expand imaging capabilities



www.amradmedical.com

SUMMIT INDUSTRIES, LLC
7555 N. Caldwell Avenue, Niles, Illinois 60714
(p) 773-353-4000 | (p) 1-800-729-9729

Certified System
SAI GLOBAL
ISO 13485
Quality



LITR044-1020