

BENNETT™

X-RAY SYSTEMS



DC-3 System

***Dedicated to Quality, Innovation and
the Chiropractic Profession***



DC-3 Film-based Radiographic System

Technology

The single most important consideration when investing in a new x-ray system is High Frequency technology. It has revolutionized chiropractic radiography. More powerful than conventional single phase or three phase generators, High Frequency technology has become the industry standard. High Frequency technology delivers the benefits of reduced patient exposure, maximum diagnostic image quality and extended tube life.

Reliability is Essential

Utilizing state of the art technology, the Bennett HF-30 generator employs internal self diagnostics and sophisticated self-testing, capable of making vital changes during the x-ray exposure to insure the integrity of the exposure is flawless.

Mechanical structures, such as tubestands and vertical cassette structures, are made from solid steel insuring years of uncompromising operation. Bennett systems also utilize stainless steel bearings for guaranteed accurate alignment of the x-ray tube and bucky.

Performance and Operation

Offering fast 2-point selection of only the kVp and mAs factors, the HF-30 generator will automatically select the highest mA station and the shortest exposure time – resulting in less image blur caused by patient motion. Anatomical Programming on the operator console (optional on the HF-30) offers consistent, effortless images by fast, automatic technique selection in three easy steps.

Digital Ready – Right For Today and Tomorrow

All of our radiographic systems are compatible with most digital imaging systems. Our IMAGER Upgrade may be added to your system any time after purchase, converting it to a DICOM-compliant digital imaging system.

Built-in Quality Protects Your Investment

Our rigorous manufacturing process and adherence to the highest level quality assurance measures guarantees that every product will provide years of trouble-free operations. As a testimonial to our commitment to quality, Bennett proudly backs every chiropractic film-based radiographic system with a 7-year parts warranty, the best warranty in the industry.



**Optional Traditional
Wallstand**

BENNETTTM

2901 W. Lawrence Avenue • Chicago, IL 60625
800-641-4107 • 773-588-2444 • Fax 773-588-3424
www.bennettxray.net