

# Radiographic Floating Top Tables



A4000 Elevator table

**AMERICOMP**

# All the Features You Need at a Price You Can Afford


Since 1987, the quality standards of Americomp have been unsurpassed. Americomp radiographic components remain the affordable choice in today's marketplace, offering the same long-term dependable performance that has become synonymous with the Americomp name. Our manufacturer's parts warranty secures your investment today and guarantees your investment value tomorrow.

## Float Top Tables

The A4000 elevating float top table provides an elevation range of 22 1/2" to 34" (57 cm to 86 cm) for effortless patient transfer. It is very well suited for immobile, pediatric or geriatric patients. The smooth, quiet lift mechanism and high-capacity motor safely handles patients weighing up to 400 lbs. (182 kg).

The A5480 float top table is 31" (79 cm) in height with a welded steel pedestal base that facilitates patients weighing up to 400 lbs. (182 kg).

Both tables feature a generous table top workspace, allowing for comfortable and easy patient positioning. Precision stainless steel ball bearings float the table top in four directions. The heavy-duty table top locks provide a stable surface for fast and safe transfer of patients even during power outages. The recessed foot controls within the table base permit hands-free operation, while restricting inadvertent patient release.



The large, flat working surface of the table provides images without distortion.



A5480 Four way Table

## High Performance Generators

Complement any Americomp radiographic table with a high-quality, state-of-the-art generator. Contact your Americomp dealer for details.

## Digital-ready

Americomp radiographic components are is completely compatible with most digital imaging systems. Digital imaging technology may be added at any time.



Quality  
ISO 9001  
SAI GLOBAL

## Americomp Inc.

Chicago, Illinois, USA

1-800-458-1782

Fax 773-588-3424

[www.americompus.net](http://www.americompus.net)

LITEGO-0713