


* These drawings have been provided to use as a conceptual placement of equipment. This is not considered a valid room drawing.



EQUIPMENT LEGEND					
FUTURE					
EXISTING					
INSTALLED BY CUSTOMER/CONTRACTOR					
FURNISHED BY CUSTOMER/CONTRACTOR					
INSTALLED BY SUMMIT					
FURNISHED BY SUMMIT					
ITEM	DESCRIPTION DIM: W X D X H	WT. (LBS)	BTU'S /HR	DETAIL	
①	Summit SHFR Generator Cabinet DIM: 17.5" x 16.9" x 22.2"	143	1025		• • • • •
②	Summit SHFR Generator Console DIM: 16.9" x 11.4" x 1.9"(w/Pedestal 36.6")	13.2	N/A		• • • • •
③	Wall Mounted DR Workstation PC, Monitor, and Keyboard	N/A	N/A		• • • • •
④	Summit FRS Digital Stand (Right Load) DIM: 96" x 64" x 85"	N/A	1706		• • • • •
⑤	Lead Window - size to be determined by customer	N/A	N/A		• • • • •
⑥	Summit Mobile Table DIM: 79" x 29" x 28"	101	N/A		• • • • •

Minimum Ceiling Height - 8'0"



Dwg No. S2103-05-B
Please Direct all Questions to: Susan Pryor Welter
susanw@summitindustries.net
678-993-5810

Typical Plans for:
Summit FRS System
Minimum Room Size - with Table

Page
A1
Date
3/30/21

Drawn by: Cory Kimbrough