



S109 Wallstand

Wall Stand (17 in.)

- Fully counterbalanced
- 17 in. Grid Cabinet with front panel travels from floor to 76 in.
- Grid: 10:1 ratio, 103 LPI, specify focal range
- Specify left or right load
- Heavy-duty Deluxe cassette tray 06338
- Electric locks
- Operates on 24VDC



Comfortable Chin Rest



E-Z Lock Release

TRAVEL AND ALIGNMENT CONSIDERATIONS

Wall to rear of column	1/8" (flush) or 1-1/2" (offset)
Wall to film plane	10" (flush) or 11-3/8" (offset)
Center of film plane vertical travel	72-3/4" maximum to 9-3/4" minimum

MOUNTING SPECIFICATIONS

Minimum ceiling height	85.5" from finished floor
Wall bracket mounting hole centers	Four total, 84" from finished floor
Four wall anchors, each with a minimum pull-out force of 250 lbs. Two floor anchors, with a minimum pullout force of 500 lbs.	See reference sheet for typical holding forces of various devices in common construction materials.

ELECTRICAL REQUIREMENT FOR ELECTRIC LOCK VERSION

24 VDC, 2 Amp, from D164 power supply (included) or applicable x-ray generator.