



S210 Non-Elevating Table

FOUR-WAY FLOAT TOP TABLE

This table has a 400 pound load rating. It comes with treadle operated electric locks and a flat, floating top 84 x 30 inches in size. A 76-inch top is available as an option. The table top can only be mounted by sliding it on from the end. The table requires 115 or 230 VAC input power and contains an integrated power supply for 24 VDC locks. The table base is open for easy electrical access and equipment storage if desired.

TRAVEL AND ALIGNMENT

Table top travel	Transverse: ± 5 " from center Longitudinal: ± 14 " from center
Table top height	31.5 inches
Patient load rating	400 pounds
Rear of upper table base to film center	14.0 ± 0.2 inches

ELECTRICAL REQUIREMENTS

115VAC @ 1.25A/230VAC @ 0.6A, 50/60 Hz
A 20 foot line cord with IEC plug is provided.

MOUNTING SPECIFICATIONS

The table bolts to the floor at four places. Recommended holding force for each bolt (not provided) is 500 pounds.



SHIPPING INFORMATION

CONTENTS	DIMENSIONS	WEIGHT
Table top (Standard 84")	4 x 34.5 x 87.5 inches	80 lbs. (est.)
Table base	38 x 36 x 56 inches	430 lbs. (est.)

Specifications subject to change without notice.

FOUR-WAY TABLE DRAWING

S210 with 84" Table Top



SUMMIT INDUSTRIES, LLC
7555 N. Caldwell Avenue, Niles, Illinois 60714
(p) 773-353-4000 | (p) 1-800-729-9729