



CS Classic

WALL MOUNTED CHEST STAND SYSTEM



The TXR International **CS Classic Chest Stand System** accommodates a full range of upright exams and is designed for facilities that specialize in pulmonology and occupational medicine. The CS Classic System comes with an industry leading 5-year parts warranty.

The extensive vertical travel accommodates a full range of upright exams.

Specifications

Tube Stand (S114)

- Floor-to-wall or floor-to-ceiling mount
- Angulation dial for easy positioning
- 10" transverse movement, non-rotational tube stand
- 360° tube angulation with detents at zero and +/-90°
- 13.25 in. minimum (focal spot to floor), 73.5 in. maximum vertical travel
- Floor-to-ceiling mount for 88 in. – 121 in. ceiling height

CS Classic Wall Stand (S109)

- Vertical travel of detector center: 10" – 72" above finished floor
- Fully counter-balanced
- Electric locks
- Accepts all size detectors

Generator

- 30 - 50kW generator options

Options

- Pair with optional mobile table to expand imaging capabilities

Specifications are subject to change without notice




DR UPGRADEABLE
ASK ABOUT OUR DR UPGRADE SOLUTIONS




Chin rest



EZ lock release



is a Trusted Brand of Summit Industries



Since 1984

SUMMIT INDUSTRIES, LLC
7555 N. Caldwell Avenue, Niles, Illinois 60714
(p) 773-353-4000 | (p) 1-800-729-9729
www.summitindustries.net



LITX001-0323
CS Classic