



FWFC *Classic*

FLOOR-TO-WALL / CEILING RADIOGRAPHIC SYSTEM

Floor-to-Wall / Ceiling Tubestand

- The full-featured FWFC can be configured to meet the needs of any department
- Cross-table exams are easily performed using the $\pm 180^\circ$ rotational and ± 5 in. transverse positioning capabilities
- Available with either 8 ft. or optional 10 ft. rail, the FWFC can be used to perform weight bearing studies



Product Highlights



Compact design for facilities with space constraints



Economical and versatile to meet the budgetary requirements of the facility



Customized and flexible configurations



Industry Leading Warranty

Four-way Float Top Table

- Flat, four-way, floating top permits easy transfer and positioning of patients
- Designed for all types of applications, available in 400 lb patient capacity
- Recessed foot-treadle allows for safe patient transfer



Wallstand

- Designed for ease of use and simplicity, the wallstand lowers for weight bearing exams like standing knees
- Ergonomically designed handle allows for easy positioning of the wall receptor
- The wallstand can be used with CR and DR panels



*Specifications subject to change without notice.



www.amradmedical.com

SUMMIT INDUSTRIES, LLC
7555 N. Caldwell Avenue, Niles, Illinois 60714
(p) 773-353-4000 | (p) 1-800-729-9729



LITR033-0523