

Attention: Sales Staff

AMERICOMP Dealer Notice

March 19, 2024

More good news ! Skinny Chiropractic Systems are back.

Thin is in! Due to dealer requests we are continuing this program to provide what you need: basic systems with options allowing your customer to only pay for the features they need.

See the following pages for details.

For a limited time only, we are pleased to provide the **industry's best value** on three **Skinny AC1se** systems:



- **20kW HF system : \$ 9,937 net !**
- **30kW HF system : \$ 10,087 net !**
- **40kW HF system : (1 phase or 3 phase) \$ 11,480 net**

To take advantage of this offer, orders must be placed and shipped by **May 31, 2024.**

Other promotions do not apply.

There is additional money to be made on options, read on for details.

Hurry !don't miss out.



SHOP AMERICOMP AND SAVE MONEY!

Call us Toll Free 1-800-458-1782 or 773-353-4000

AMERICOMP

QUOTED TO: All Americomp Dealers		DATE: March 19, 2024
QTY.	PRODUCT QUOTED	NET PRICE
1	<p>Skinny AC1se System consisting of:</p> <p>Select a Generator:</p> <ul style="list-style-type: none"> - Spectra 20HF, 20 kW Generator <ul style="list-style-type: none"> > without 13 ft. line cable (see options to add) > 15 ft. console cables (see options for 50 ft.) > service manual on CD (see options to add printed version) > without UL (see options to add) > Operates on 1 phase, 208 VAC, 70 Amp or 240 VAC 60 Amp electrical service - Spectra 30HF, 30 kW Generator <ul style="list-style-type: none"> > without 13 ft. line cable (see options to add) > 15 ft. console cables (see options for 50 ft.) > service manual on CD (see options to add printed version) > without UL (see options to add) > Operates on 1 phase 240 VAC, 100 Amp electrical service (see next page for options for other electrical services) - Spectra 40HF AP, 40 kW HF Generator <ul style="list-style-type: none"> > without 13 ft. line cable (see options to add) > 15 ft. console cables (see options for 50 ft.) > service manual on CD (see options to add printed version) > without UL (see options to add) > Operates on 1 phase, 240 VAC, 150 Amp or 1 phase 208 VAC, 175 Amp or 3 phase, 240 VAC, 100 Amp electrical service. <p>A4570 – Tubestand</p> <ul style="list-style-type: none"> - Floor-to-wall or ceiling (specify at time of order), counterbalanced column - 5-foot tracks, manual lock - Angulating tubemount with manual lock - Specify tube arm length: 10", 14", 17", 19" or 23" (tube focal spot to front face of tubestand) - Angulation dial and operator handles <p>A120 – Tube</p> <ul style="list-style-type: none"> - 140,000 H.U., 1.0 – 2.0 mm focal spot (specify shoot left or right) <p>C600-15 – Cables</p> <ul style="list-style-type: none"> - 15' high voltage cables <p>07877-125 – Collimator</p> <ul style="list-style-type: none"> - Certified, manual collimator with LED, dual laser crosshairs and swivel <p>Does not include:</p> <p>Wallstand, grid cabinet, grid or tray – feel free to add a new one from the options, or your own used one, or add a digital product.</p> <p>Parts Warranty: 1-year</p> <p>Terms:</p> <p>Prices are FOB Niles, IL and do not include shipping or taxes. Payment terms: Net 30 days with approved credit. Prices are valid for orders placed and shipped by May 31, 2024.</p>	<p>\$9,937 each</p> <p>\$ 10,087 each</p> <p>\$ 11,480 each</p>
	See page 3 for options.	

QUOTED TO: All Americomp Dealers		DATE: March 19, 2024
QTY.	PRODUCT QUOTED	NET PRICE
	Skinny AC1se Options: Generator Options: <ul style="list-style-type: none"> ➤ For operation from 1 phase, 208-277 VAC (30HF only) 08168-001 add: \$ 275 ➤ For operation from 3 phase, 240 VAC, 60 Amp or 3 phase, 480 VAC, 35 Amp electrical service (30 HF only and includes AP & 500 mA) add: \$ 1,097 ➤ For operation from 3 phase, 208-240 VAC electrical service (for 40 HF & 50 HF only) 08168-002 add: \$ 817 ➤ For operation from 3 phase, 380-480 VAC electrical service (for 40 HF & 50 HF only) 08168-003 add: \$ 817 ➤ Anatomical programming upgrade (20HF & 30HF only) add: \$ 487 ➤ UL label add: \$ 58 ➤ Printed service manual in binder add: \$ 58 ➤ 13 ft. line cable (from power module to electrical disconnect) add: \$ 81 ➤ 500 mA upgrade (30HF only) add: \$ 395 ➤ 50 ft. console cable upgrade add: \$ 58 ➤ Wall mount for operator console L311 add: \$ 113 ➤ Pedestal mount for operator console L310 add: \$ 375 ➤ DR software integration add: Call for compatibility Tubestand Options: <ul style="list-style-type: none"> ➤ Electric locks upgrade add: \$ 213 ➤ 8 ft. rails replace 5 ft. rails add: \$ 126 ➤ 10 ft. rails replace 5 ft. rails add: \$ 244 ➤ Fixed tubestand A3284 (vertical tube travel only) replaces variable tubestand A4570 (vertical and horizontal tube travel) deduct: (\$ 84) ➤ Transverse Tubestand S114, electric locks, 8 ft. rails replaces A4570 Tubestand (we recommend adding 20 ft. H.V. cables) add: \$ 1,429 Wallstand Options: <ul style="list-style-type: none"> ➤ Wallstand A1831, manual lock, 17" grid cabinet, front panel, 8:1 103 LPI grid, regular-duty tray add: \$ 1,042 ➤ Electric lock upgrade for above Wallstand add: \$ 193 ➤ Heavy Duty Cassette Tray 03092 upgrade for above Wallstand add: \$ 188 ➤ Wallstand only A3546, manual lock, <u>without</u> grid cabinet, grid or tray add: \$ 604 ➤ Wallstand A1831, electric lock, w/DR panel holder 06778 add: \$ 1,234 ➤ Tilt Wallstand S102, electric lock, 17 " grid cabinet, front panel, 10:1 103 LPI grid, Heavy-duty Tray 03092. We recommend adding 10 ft. tubestand rails and 25 ft. H.V. cables to take extremity views (cabinet horizontal) add: \$ 3,499 Option: DR panel holder 06778 replaces grid cabinet, grid & tray add: n/c Other Options: <ul style="list-style-type: none"> ➤ Mobile table 24" x 76" w/grid cabinet, 8:1 103 LPI grid, regular duty tray add: \$ 1,688 ➤ 200,000 H.U., 0.6-1.5 mm FS tube upgrade add: \$ 435 ➤ 10:1 103 LPI grid upgrade add: \$ 76 ➤ Additional 5 ft. of H.V. cables add: \$ 58 Warranty Options: <ul style="list-style-type: none"> ➤ Extend parts warranty from 1-year to 5-years, x-ray tube pro-rated add: \$ 487 total 	