


TS Primo DR

DIGITAL IMAGING PACKAGE




Digital Radiography Made Easy



Single-Screen Setup

Navigate easily with our all-in-one screen setup menu. Less clicks to find the information you need.



One-Screen Registration

Convenient one-screen registration makes entering emergency cases quick and easy. Worklist feature is compatible with all HIS / RIS systems.



Easy Acquisition

Populate exams from a wide selection of views. Add studies, change and delete with efficiency. Suspend studies for later completion or send to PACS.

Ease of Use, Peace of Mind



Putting Innovation in Your Toolbox

Having the right tools in your tool box helps complete jobs quickly and efficiently. AMRAD Medical ATS Primo DR has all of the tools from measuring and angle tools to crop and mask all on one screen, streamlining the acquisition process.

Smart UI Touch

The Power of Touch


24" Touchscreen monitor puts the power of DR at your fingertips. Large screen delivers bright visuals in high definition. Sensitive touchscreen with faster selections for improved productivity.

Fast and intuitive post processing



Detector Specifications

1417



1717



Detector Technology	Amorphous Silicon	Amorphous Silicon
Scintillator	CsI	CsI
Active Area	14" x 17"	17" x 17"
Number of Pixels	3500 x 4300	4267 x 4267
Pixel Pitch	100 um	100 um
AD Conversion	16 bit	16 bit
Trigger Mode	AED (Optional) / Software	AED (Optional) / Software
Data Interface	DICOM connections using Gigabit wired or IEEE 802.11n (2.4 and 5GHz) wireless Ethernet	DICOM connections using Gigabit wired or IEEE 802.11n (2.4 and 5GHz) wireless Ethernet
Dimension	384mm x 460mm x 15mm	460mm x 460mm x 15mm
Weight	6.6 lbs.	7.5 lbs.
Water Proof	IP56	IP56
Drop Height	100cm @3mm PVC	100cm @3mm PVC
Static Loading	661 lbs. Uniformly	661 lbs. Uniformly
Battery Capacity	8.5 h	8.5 h
Battery Charging	4 h in Charger	4 h in Charger

Workstation and Monitor Specifications:

- L-P360 (based on availability, model + specs may vary)
- 14th Generation Intel® Core™ i3-14100 Processor (P-cores 3.50 GHz up to 4.70 GHz)
- 16 GB DDR5-4400MHz
- Intel UHD Graphics
- 1 TB SSD M.2 2280 PCIe Gen4
- C7-All SATA Hard Drives, RAID1 FOR 2 Hard Drives
- OS installed on a 256 GB SSD M.2 2280 PCIe Gen4
- Windows 10 Enterprise LTSC 64 Bit OS
- 16X DVD+/-RW
- 5-Year Parts warranty
- 24-Inch Full HD 1080p (1920x1080) 16:9 LED Touch Screen Monitor with HDMI 1.4, DisplayPort, USB 3.0, and VGA Inputs, robust, scratch resistant 10-point multi-touch full HD 1080p touchscreen, 50M:1 Dynamic Contrast Ratio

All specifications are subject to change without notice.



SUMMIT INDUSTRIES
A NETWORK OF SERVICE & SUPPORT

Summit Industries' successful service network surrounds you with a dynamic team of exceptional experts. Your organization and your patients will benefit from the experience and best practices that only Summit Industries and their Authorized Dealer Network can offer.

AMRAD  **Medical**
by Summit Industries

SUMMIT INDUSTRIES, LLC
7555 N. Caldwell Avenue, Niles, Illinois 60714
(p) 773-353-4000 | (p) 1-800-729-9729

LITR048-1124



www.amradmedical.com