



# CS Elite

WALL MOUNTED CHEST STAND SYSTEM



The AMRAD Medical CS Elite system is designed for facilities that specialize in Pulmonology and Occupational Medicine. Advanced features like ergonomically correct release handle, electric locks and remote collimator light make positioning fast and efficient.

- 30-50kW generators
- Accommodates a full range of upright exams
- Pair with an optional mobile table to expand imaging capabilities
- DR upgradeable - ask about our digital upgrade solutions


## Specifications

### Floor-to-Wall or Floor-to-Ceiling, Non-Rotating with Transverse Travel Tubestand

- Fully counterbalanced
- Platform tube mount
- 24 in. angulating tube arm
- +/- 5 in. transverse tube travel
- Angulation dial and operator handgrips
- Vertical, horizontal, transverse and tube arm angulation electric locks
- 8 ft. rails (provides 6.5 ft. travel)
- Floor-to-ceiling mount for 88 in. – 121 in. ceiling height
- 13.25 in. minimum (focal spot to floor), 73.5 in. maximum vertical travel
- Operates on 24VDC

### Wallstand

- Versatile center-mounted or flag mounted space saving design
- Fully counterbalanced
- 17 in. Grid Cabinet with front panel
- Grid: 10:1 ratio, 103 LPI, specify focal range
- Heavy-duty Deluxe cassette tray 06338
- Electric locks
- PA Grips (2 with center mount style, 1 with flag mount style)
- Operates on 24VDC
- Remote Collimator Light



**SUMMIT INDUSTRIES**  
A NETWORK OF SERVICE & SUPPORT

Summit Industries' successful service network surrounds you with a dynamic team of exceptional experts. Your organization and your patients will benefit from the experience and best practices that only Summit Industries and their Authorized Dealer Network can offer.



**SUMMIT INDUSTRIES, LLC**  
7555 N. Caldwell Avenue, Niles, Illinois 60714  
(p) 773-353-4000 | (p) 1-800-729-9729



www.amradmedical.com

LITR045-0225