



SRS Pocket SMALL ROOM SOLUTION
FLOOR-TO-WALL MOUNTED TUBESTAND



Big Room Capability That Fits in a Small Pocket of Space!!

Small room size? Need more flexibility within that small space?

The AMRAD Medical SRS Pocket is the solution! The SRS Pocket is a compact system that was specifically designed to fill the needs of facilities with limited space. This system can fit into a space as small as 9' x 9', while providing the functionality of a full conventional X-ray system.

The SRS Pocket has a multi-functional interlocking tube stand and wall stand with a dual rail design that makes positioning easy and efficient. When the need arises, the tube stand and wall stand can be separated to further improve the functionality of such a small space.



Perform a full range of exams including:

- ✓ Upright chest
- ✓ Extremities
- ✓ Weight bearing studies
- ✓ Table exams (with optional mobile table)

Peace of Mind

The SRS Pocket is backed by our industry leading 5-year parts warranty.



Seeing is Believing

The AMRAD SRS Pocket is designed with our “EZ-Read” angulation dial, taking the guesswork out of getting that exact angle. Specially designed detents for 40” to table and 40” and 72” to wall make positioning quick and efficient. Our “SSD” system status display lets you know when you are in the detents and ready to perform your exam, improving throughput and productivity.



The easy-read backlit blue angulation dial makes identifying and adjusting angles fast and accurate.



Tilt the Odds in Your Favor

A tilting grid cabinet with detents every 15° that is available as a fixed DR cabinet, or with a removable tray for wireless panels, adds even more flexibility. Easy one-handed rotation of the grid cabinet allows precise positioning for exams such as cervical spines or sunrise knees.



SRS Pocket Specifications

Wall Stand

WALL-MOUNTED WALL STAND (17" X 17")

- Tilting 17-inch grid cabinet with front panel that can be positioned from 0° (horizontal) to 105° (15° past vertical) with 15° locking increments
- 10:1 103 LP Grid
- Fixed DR grid cabinet (see receptor tray in options)
- Electric locks: horizontal and vertical
- Remote collimator light switch
- Fully counterbalanced
- Operates on 24VDC powered by Nexus box

Generators

- Choose from our 40 kW - 80 kW HF generator options



Tube Stand

FLOOR-TO-WALL RAIL-MOUNTED TUBE STAND

- Platform tube mount
- Fixed length angulating tube arm
- Angulation dial and operator hand grips
- Vertical travel: 13.5" minimum to 73" maximum (focal spot to floor)
- Electric locks: vertical, horizontal, and tube arm angulation
- Fully counterbalanced
- Operates on 24VDC powered by Nexus box

Mobile Table (Optional)

- 76 in. x 24 in. flat, seamless, table top (56 in. x 24 in. optional)
- 20 in., 24 in. or 28 in. table height (specify at time of ordering)
- Hospital-grade swivel casters with two locks
- Electrical power not required
- 300 lb. rated patient load



SUMMIT INDUSTRIES
A NETWORK OF SERVICE & SUPPORT

Summit Industries' successful service network surrounds you with a dynamic team of exceptional experts. Your organization and your patients will benefit from the experience and best practices that only Summit Industries and their Authorized Dealer Network can offer.



LITR047-0225

SUMMIT INDUSTRIES, LLC

7555 N. Caldwell Avenue, Niles, Illinois 60714

(p) 773-353-4000 | (p) 1-800-729-9729



www.amradmedical.com