

Proven Performance in Chiropractic Imaging

The Choice for Simplicity, Reliability and Value




AC1 17" Chiropractic System with options

AMERICOMP

All the Features You Need at a Price You Can Afford

Since 1987, the quality standards of Americomp have been unsurpassed. Americomp systems remain the affordable choice in today's marketplace, offering the same long term dependable performance that has become synonymous with the Americomp name.



AC1 Chiropractic System Features

Ideal for an existing practice or start-up practitioner, this reliable and versatile system is an outstanding value.

The AC1 system accommodates up to a 17" (43 x 43 cm) cassette or digital detector in either direction, allowing a broad range of procedures.

The system features a traditional floor-wall or ceiling tubestand, and provides film focal distances of 40" to 72". The tubestand is constructed of heavy-duty steel, insuring many years of reliable service. Stainless steel ball bearings allow effortless positioning of the x-ray tube and precise alignment to the wallstand over the entire range of travel. The tubestand angulation dial is an easy-to-read, accurate indicator of tube angle. Sturdy locks provide secure placement of the x-ray tube. Optional electric locks allow quick tube positioning.

The wallstand utilizes a precision bearing system that provides effortless vertical travel of the grid and cassette tray. The full travel wallstand allows a wide range of exams including standing knee and ankle procedures. The wallstand features a positive manual lock or optional electric lock.

Optional S102 tilting wallstand allows for upper cervical and easy extremity imaging.



AC1 with optional S102 tilting wallstand and extended rails for extremity imaging

High Performance Generators

Complement any Americomp radiographic system with a high-quality, state-of-the-art generator. Contact your Americomp dealer for details.

Digital-ready Systems

Americomp radiographic equipment is completely compatible with most digital imaging systems. Digital imaging technology may be added at any time.

Americomp
7555 N. Caldwell Ave.
Niles, Illinois 60714
1-800-458-1782
www.americompus.net